| | | |
|----|---------------------------|------------------------------------|
| EN | Instructions for use | PTCA Balloon Dilatation Catheter |
| FR | Mode d'emploi | Cathéter de dilatation pour ACTP |
| DE | Gebrauchsanweisung | Ballondilatationskatheter (PTCA) |
| IT | Istruzioni per l'uso | Catetere dilatatore a PTCA |
| NL | Gebruiksaanwijzing | PTCA Ballon Dilatatie Katheter |
| ES | Instrucciones de Uso | Catéter de dilatación (ACTP) |
| CS | Návod k použití | Balónkový dilatační katetr (PTCA) |
| EL | Οδηγίες χρήσης | Καθετήρας μπαλόνι διαστολής (PTCA) |
| PL | Instrukcja stosowania | Balonowy cewnik dylatacyjny (PTCA) |
| PT | Instruções de Uso | Cateter Balão de Dilatação (PTCA) |
| TR | Kullanım talimatlarına | PTCA Balon Dilatasyon Kateteri |
| UK | Інструкція з експлуатації | Катетер для балонної дилатації |
| ZH | 中文 使用说明 | PTCA 球囊扩张导管 |

EVEREST

PTCA CATHETER

SUMMIT

CTO CATHETER

FORCE NC

NON-COMPLIANT

w
welling®

Percutaneous Transluminal Coronary Angioplasty (PTCA) Balloon Dilatation Catheter.

Description

The distal section is dual lumen, one lumen permits the inflation of the balloon, the second lumen permits the use of a 0.014" (0.36 mm) guide wire, to facilitate the advance of the catheter to and through the stenosis. The proximal section of the catheter is a single lumen, stainless steel hypotube with a single luer port for connection to the inflation/deflation device.

The balloon includes one or two radiopaque markers to aid in the positioning of the balloon within the stenosis. Marks located on the proximal portion of the catheter shaft indicates the exit of the balloon catheter tip out of the guide catheter (one at 90cm and one at 100cm).

Intended Purpose / Intended Use

The PTCA balloon dilatation catheter is intended for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial reperfusion in patients requiring a percutaneous coronary interventional (PCI) procedure in a blood vessel. The PTCA balloon dilatation catheter is intended for transient use.

Indications

- Indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial reperfusion.

Contra-Indications

The use of the balloon dilatation catheter is contraindicated for:

- Unprotected left main coronary artery
- Coronary artery spasm in the absence of a significant stenosis

Patient Target Groups

The PTCA catheter is intended to be used in patients suffering from stenotic lesions in the coronary vascular system who require a PCI.

Intended User(s)

The target users are medical professionals (e.g. Interventional Cardiologists, cardiovascular technicians) who perform and assist in catheterisation procedures in a clinical laboratory setting, during percutaneous intervention (e.g. PTCA, coronary stent placement).

Use of the product is restricted to healthcare professionals in an interventional suite with angiography equipment, with personnel with relevant and adequate training and who are familiar with the possible/conceivable complications.

Clinical Benefits

PTCA balloon catheter characteristics for procedure technical success following successful guidewire crossing of the target lesion, with the benefit of optimizing lesion preparation and/or postdilatation for an overall successful PCI procedure for improving myocardial reperfusion.

PTCA balloon catheter characteristics for procedure technical success of lesion preparation and/or post-dilatation following successful guidewire crossing of the target lesion (as part of an overall PCI procedure), thereby contributing to providing relief of the patient's coronary ischemia symptoms.

Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse as this can potentially result in compromised device performance and increased risk of inappropriate reesterilization and cross contamination.
- To reduce the potential vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are acceptable candidates for coronary artery bypass graft surgery requires careful consideration. This includes possible hemodynamic support during PTCA as treatment of this patient population carries special risk.
- When the balloon dilatation catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic guidance. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon Pressure: do not exceed the rated burst pressure (RBP) indicated on the packaging. The rated burst pressure (RBP) is based on results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure (RBP).
- Use of a pressure-monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary bypass surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to "Use Before" ("expiry") date specified on the package.
- Do not use if packaging is damaged or unintentionally opened prior to use
- Treatment of vessels with moderate or heavy lesion calcification is associated with decreased success rates up to 60–85% and increases the risk of acute closure, vessel trauma, balloon burst, balloon entrapment and associated complications.

Precautions

- Prior to angioplasty procedure, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it will be used.
- Only physicians trained in the performance of percutaneous transluminal coronary angioplasty should use the catheter.
- During the procedure, appropriate anticoagulant and coronary vasodilator therapy must be provided to the patient as needed. Anticoagulant therapy should be continued for a period of time to be determined by the physician after the procedure.

Adverse effects

Possible complications include, but are not limited to the following by performing of a Percutaneous Transluminal Coronary Angioplasty procedure:

- Death.
- Acute myocardial infarction.
- Total occlusion of the coronary artery or bypass graft.
- Coronary vessel dissection, perforation, rupture or injury.
- Restenosis of the dilated vessel.
- Haemorrhage or haematoma.
- Unstable angina.
- Arrhythmias, including ventricular fibrillation.
- Drug reactions, allergic reaction to contrast medium.
- Hypo/hypertension.
- Infection.
- Coronary artery spasm.
- Arteriovenous fistula.
- Embolism.

Inspection prior to use

Prior to use, carefully examine the catheter packaging (pouch & box) for any damage. Prior to use of the balloon dilatation catheter in PTCA, all equipment to be used for the procedure, including the dilatation catheter should be carefully examined for defects. Examine the balloon dilatation catheter for bends, kinks or other damage. Do not use any defective equipment. Prepare equipment to be used following manufacturer's instruction or standard procedure.

Preparation of the Balloon Dilatation catheter

Complete the following steps to prepare the catheter:

- Carefully remove the balloon dilatation catheter from the pouch and place in/onto a sterile field.
- Remove balloon dilatation catheter from hoop.
- Gently remove the balloon protector together with the stylet from the balloon dilatation catheter. Prepare an inflation device with the recommended contrast medium according manufacturer's instruction.

To evacuate air from the balloon segment, the following procedures should be followed.

- Fill a 10cc syringe with approximately 2cc 50/50–mixture contrast/saline.
- Connect the 10cc syringe to the balloon dilatation catheter, hold the nozzle of syringe downwards and pull the plunger to create a negative pressure (hold this for at least 15 seconds).
- Slowly release plunger to neutral (with nozzle of syringe pointed downwards), allowing contrast to fill the shaft of the balloon dilatation catheter.
- Disconnect syringe from balloon dilatation catheter and remove air out of syringe.
- Repeat steps 2 to 4 until no more bubbles appear in the syringe.
- Disconnect syringe from the balloon dilatation catheter and connect the inflation device, with a fluid–fluid connection, to the balloon dilatation catheter.
- Pull the plunger of the inflation device to create a negative pressure.

Warnings

- An incorrectly prepared balloon dilatation catheter may prolong the deflation time.
- Loss of vacuum or a continuous stream of air entering the syringe upon aspiration indicates the presence of a leak in the system (inspect the system for integrity).
- Do not torque the balloon dilatation catheter shaft.

Insertion technique

- First advance the guiding catheter and connect to Y–connector (minimal 0.056" / 1.42 mm ID) and cross the stenosis with a 0.014"/ 0.36 mm (or smaller) guide wire according to the manufacturer's instructions.
- Wipe the extending guide wire with a sponge or gauze to remove residuals of blood or contrast.
- Load the guide wire into the distal tip of the balloon dilatation catheter guide wire lumen and ensure that the guide wire exits at the transition point 25 cm proximal of the distal tip.
- Advance the balloon dilatation catheter over the guide wire until it approaches the Y–connector.
- Loosen the valve of the Y–connector and insert the balloon dilatation catheter for minimal 30 cm while maintaining the guide wire position.
- Close the Y–connector to create a seal around the balloon dilatation catheter but ensure that movement is not restricted.
- Advance the balloon to the stenosis under fluoroscopy.
- Inflate the balloon while monitoring the inflation pressure, the balloon position, the ECG, blood pressure and other vital conditions of the patient.
- After dilatation of the stenosis deflate the balloon, check deflation under fluoroscopy.
- Pull the balloon dilatation catheter back, into guiding catheter, while maintaining the guide wire position.
- Perform a control angiography.
- Continue the procedure according to the accepted PTCA protocol.

Catheter exchange procedure

- Loosen the haemostatic valve.
- Hold the guide wire and haemostatic valve in one hand, while holding the balloon shaft in the other hand.
- Maintain guide wire position in the coronary artery by holding the wire stationary and pull the balloon dilatation catheter out of the guiding catheter while monitoring wire position under fluoroscopy.
- Withdraw the deflated balloon dilatation catheter until the transition point in the guide wire lumen is reached. Carefully pull the flexible, distal portion of the balloon catheter out of the haemostatic valve while maintaining the guide wire position across the lesion.
- Prepare the next balloon dilatation catheter to be used, as previously described.
- Backload the new balloon dilatation catheter onto the guide wire as previously described and continue the procedure accordingly.

Warning

After use, the equipment used in angioplasty, may be a potential biohazard. Handle carefully and dispose of in accordance with hospital policy and appropriate applicable federal regulations.

Notice to the user

In case of any serious incident that has occurred in relation to the device, this should be reported to the manufacturer and the competent authority of the Member State in which the user and/ or patient is established.

Link to the Summary of Safety and Clinical Performance (SSCP)

The Summary of the Safety and Clinical Performance (SSCP) is available in the European database on medical devices (Eudamed).

The PTCA catheters have Basic UDI-DI 872063420BMD–PTCAJ9.

<https://ec.europa.eu/tools/eudamed>

Warranty

Blue Medical Devices B.V., (BLUE), warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excluded all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties or merchantability of fitness. Handling storage, cleaning and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BLUE control directly affect the device and the results obtained from its use.

BLUE's obligation under this warranty is limited to the repair or replacement of this device and BLUE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. BLUE neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. BLUE assumes no liability with respect to devices reused, reprocessed or re–sterilized and makes no warranties expressed or implied, included but not limited to merchantability or fitness for intended use, with respect to such device.

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© Blue Medical Devices B.V.
Blue Medical Devices B.V. is an