





- Buffer solutions
- Cleaning agent
- Probe stand
- Probe extension cord
- Dust caps
- Power supply
- Batteries

Contact:

Sentron Europe BV,
member of the Welling
Group

Kamerlingh-Onnesstraat 5
9351 VD Leek
The Netherlands

sentron@welling.com
tel +31 50 501 38 00

Item	Description	Order code	Details
Buffer solutions			
	Sentron buffer solutions are N.I.S.T. based, high quality reference fluids.	A116-002	pH 2.00 buffer 0.5l +/- 0.02 pH N.I.S.T @ 25°C
	The twin neck bottle provides a convenient and easy to use way of providing the right amount of buffer without waste or chance on contamination.	A116-004	pH 4.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C
		A116-007	pH 7.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C
		A116-010	pH 10.00 buffer 0.5l +/- 0.01 pH N.I.S.T @ 25°C
		A116-012	pH 12.00 buffer 0.5l +/- 0.05 pH N.I.S.T @ 25°C
Cleaning agent			
	Cleaning agent for removal of enzymes and proteins.	A118-011	1.8kg of cleaning agent.
	Alconox Terg-A-Zyme is a concentrated powder and can be prepared by adding (deionized) water. One package makes approx. 200 liters of detergent.		
Probe stand			
	The probe stand can fit probes and includes 2 adapter rings for mounting probes with various shapes and dimensions.	A100-001	Probe stand incl. 2 adapter rings.
	The probe stand has a solid weighted base and an easy adjustable arm to which the probe cable(s) can be attached.		
Probe extension cord			
	Probe extension cord of 1.5 m.	A109-010	Probe extension cord 1.5 m.
	The extension cord provides additional length for those situations where extra distance between probe and meter is required.		

Item	Description	Order code	Details
<p>Dust caps</p> 	<p>A set of dust caps for both meter and probe connector.</p> <p>The dust caps completely seal and cover the connectors against dust and moisture when the probe and/or meter is not in use.</p>	A102-002	<p>A set of two dust caps for both the meter and the probe connector.</p>
<p>Power supply</p> 	<p>Power supply including five world plug adapters and a USB power cord.</p> <p>Adapters included are: UK, US, EU, Australia, Middle- and Far East.</p>	A102-006	<p>Power adapter set <i>in:</i> 100-240 VAC 50/60Hz, 0.2A <i>Out:</i> 5VDC , 1A</p>
<p>Batteries</p> 	<p>Replacement battery set of six AA 1.2V batteries for SI meters.</p> <p>High quality - 2400mAh NiMH - rechargeable batteries which are recommended for use in SI meters.</p>	A102-011	<p>Battery set of six batteries AA 1.2V, 2400 mAh NiMH (ANSI-1.2H2 or IEC-HR6)</p>